





Copying Cerner Auto-text into Dragon Auto-text


1. Deploy the autotext that you want to copy into Dragon into a free text area.



2. Highlight the Auto-text and select **Copy**.




3. On your Dragon tool bar, click on the **flame icon** and select manage **Auto text**.






4. Click on the plus sign on the Manage Auto-text window.



5. Paste your Auto-text into the **Content** Section.





6. Name your Dragon Auto-text as desired, then select **Apply All**.

The screenshot shows the 'Manage Auto-texts' window. On the left is a list of auto-texts, with 'review of systems' selected. The main area contains the following fields:

- Name ***: review of systems
- Description**: (empty)
- Spoken form ***: insert review of systems
- Content ***:
 - Constitutional: [No unexplained weight gain or loss, fevers, chills, fatigue]
 - Eye: [No recent visual problems]
 - ENMT: [No ear pain, nasal congestion, sore throat]
 - Respiratory: [No shortness of breath, cough]
 - Cardiovascular: [No chest pain, palpitations, leg edema]
 - Gastrointestinal: [No nausea, vomiting, diarrhea]
 - Genitourinary: [No dysuria, hematuria]
 - Hema/Lymph: [Negative for bruising tendency, swollen lymph glands]
 - Endocrine: [Negative for excessive thirst or urination, heat or cold intolerance]
 - Musculoskeletal: [No back pain, joint pain, muscle pain]
 - Integumentary: [No rash, itching, abrasions]
 - Neurologic: [No history of fainting, memory loss, numbness]
 - Psychiatric: [No anxiety, depression]
 - Allergic/Immun: [No nasal allergies, itchy/red eyes, enlarged lymph nodes]

Fields marked with * are mandatory

Apply All Close

7. To use your Dragon Auto-text, go to a free text area and command your Dragon mic, "insert (the name of your Auto-Text)"